



Performance
Made Compact

The next generation of Dynamic Support Technology™

The F&P Evora™ Full mask is uniquely designed to optimise performance throughout the night.

Dynamic Support Technology combines Stability Wings with the Floating Seal to allow for freedom of movement while keeping the mask in place.

Stability



rated the Evora Full mask as stable or very stable¹

Performance



rated the Evora Full seal performance as very good or good²




Secure Comfort.

The decoupled seal accommodates for movement while the soft silicone gently wraps under the nose for a comfortable night's sleep.

Comfortable



rated the Evora Full mask as comfortable or very comfortable³



Considering patient preference

Up to 60% of PAP therapy users report issues such as poor mask fit, air leak, redness or skin irritation, and claustrophobia.⁴⁻¹⁰

When choosing a new mask, consider patient preference.



never or rarely experienced unintentional air leak into their eyes while using the Evora Full mask¹¹



were able to sleep in their preferred sleeping position while using the Evora Full mask¹²



never or rarely felt claustrophobic while using the Evora Full mask¹³



F&P Evora Full

Sizing options



Extra Small



Small to Medium



Large



Std.

Standard

XL

Extra Large

Fitting summary



1. Adjust the top strap



2. Adjust the bottom strap



3. Check the space

For further instruction, please refer to the Fitting Your Mask section of the Evora Full mask Use and Care Guide.

Product Codes

Code	Description
EVFIXA	Evora Full Mask - Extra Small
EVFIMA	Evora Full Mask - Small to Medium
EVFILA	Evora Full Mask - Large
EVFIXMLU	Evora Full Mask Fitpack - Extra Small, Small to Medium, Large

Spare Parts

Code	Description
400EVF111	Evora Full Seal - Extra Small Spare
400EVF112	Evora Full Seal - Small to Medium Spare
400EVF113	Evora Full Seal - Large Spare
400EVF121	Evora Full Headgear Spare - Standard
400EVF122	Evora Full Headgear Spare - Extra Large
400EVF141	Evora Full Headgear Clips & Frame Clip
400EVF151	Evora Full Tube Spare
400EVF161	Evora Full Frame Spare
400FPH341	F&P Swivel White

F&P Evora Full is available on the F&P myMask™ App

Instructs patients on mask use through step-by-step fitting, fine-tuning and cleaning videos, download the F&P myMask App. Scan the QR code, or download directly.



1. 40 out of 44 participants rated the Evora Full stability as either very stable or stable. Internal validation study conducted with 44 participants in New Zealand 2020. 2. 39 out of 43 participants rated the overall seal performance of the Evora Full mask as very good or good. Internal validation study conducted with 44 participants in New Zealand 2020. 3. 41 out of 44 rated the Evora Full mask as comfortable or very comfortable. Internal validation study conducted with 44 participants in New Zealand 2020. 4. Ewles et al 2018. doi:10.1183/13993003.congress-2018.pa2372. 5. Patout et al 2018. DOI: 10.1183/13993003.congress-2018.pa2378. 6. Andreu et al 2012. doi:10.1183/09031936.00013311. 7. Brown et al 2013. doi:10.1016/j.jstrokecerebrovasdis.2011.06.010. 8. Noseda et al 2009. doi:10.1016/j.sleep.2008.08.010. 9. Maso et al 2004. doi:10.1164/rccm.200312-1787OC. 10. Nussbaumer et al 2006. doi:10.1378/chest.129.3.638. 11. 41 out of 44 participants never or rarely experienced unintentional mask leak into the eyes when sleeping with the Evora Full mask. Internal validation study conducted with 44 participants in New Zealand 2020. 12. 42 out of 44 participants were able to sleep in their preferred sleeping position while using the Evora Full mask. Internal validation study conducted with 44 participants in New Zealand 2020. 13. 43 out of 44 participants never or rarely felt claustrophobic while wearing the Evora Full mask. Internal validation study conducted with 44 participants in New Zealand 2020.

Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Google Play and Google Play logo are trademarks of Google LLC. F&P, Evora, myMask, Dynamic Support Technology and VentiCool are trademarks of Fisher & Paykel Healthcare Limited. For patent information, please see www.fphcare.com/ip. REF 622069 REV A 2021-08 © 2021 Fisher & Paykel Healthcare Limited.

www.fphcare.com

Fisher & Paykel
HEALTHCARE